


TAG
STRAIGHT




TAG STRAIGHT



Platform switching
1 mm machined neck

- Self-tapping straight dental implant
- Internal hexagon connection
- Single platform for all diameters
- SAE Surface Treatment for best osseointegration process


The internal hexagon connection, with conical implant-prosthetic support, ensure greater stability reducing micro movements between the fixture and the abutment.




Ø	L. 7	L. 8,5	L. 10	L. 11,5	L. 13
	mm	mm	mm	mm	mm
3.8	-	TAG3.8X8.5ST	TAG3.8X10ST	TAG3.8X11.5ST	TAG3.8X13ST
4.25	TAG425X7ST	TAG425X8.5ST	TAG425X10ST	TAG425X11.5ST	TAG425X13ST
5.0	-	TAG5X8.5ST	TAG5X10ST	TAG5X11.5ST	TAG5X13ST

Possibility of fixing at different transgingival heights

1mm out of bone crest



























0,5mm out of bone crest



0mm out of bone crest



	Healing cap		Transfer open tray		Transfer closed tray
	H.2 - cod. TAGVG2 H.4 - cod. TAGVG4		cod. TAGTR		cod. TAGTRAC
	Analog		Straight abutment		Straight abutment
	cod. TAGANA		H.1 - cod. TAGMD1 H.2 - cod. TAGMD2 H.3 - cod. TAGMD3 H.4 - cod. TAGMD4		cod. TAGMD
	Zero abutment		Friction fit abutment		Anatomical straight abutment
	cod. TAGMD0		cod. TAGMDFF		H.1 - cod. TAGMDA1 H.2 - cod. TAGMDA2
	Tbase abutment NR		Moncone Tbase Rot.		Angled abutment 15° cod. TAGMA15
	H. 0.5 - cod. TAGTB05 H. 1 - cod. TAGTB1 H. 2 - cod. TAGTB2		H. 0.5 - cod. TAGTBR H. 1 - cod. TAGTBR1 H. 2 - cod. TAGTBR2		Angled abutment 25° cod. TAGMA25
	Anatomical angled abutment 15°		Anatomical angled abutment 25°		Titanium temporary abutment
	H.1 - cod. TAGMAA151 H.2 - cod. TAGMAA152		H.1 - cod. TAGMAA251 H.2 - cod. TAGMAA252		cod. TAGMP
	Titanium temporary abutment non engaging		Castable cylinder		Castable cylinder non engaging
	cod. TAGMPR		cod. TAGCC		cod. TAGCCR
	Base Cr/Co		Base Cr/Co non engaging		Tag lok
	cod. TAGCR		cod. TAGCRR		H.1 - cod. TAGLC1 H.2 - cod. TAGLC2 H.3 - cod. TAGLC3 H.4 - cod. TAGLC4 H.5 - cod. TAGLC5
	Ball abutment		Open low cap with o-ring		Teflon cap cod. TAGCAPT
	H.1 - cod. TAGASF1 H.2 - cod. TAGASF2 H.3 - cod. TAGASF3 H.4 - cod. TAGASF4		cod. TAGCAP		Cap basket cod. TAGCONT

	<p>Straight Mua</p> <p>H.1 - cod. TAGMUA1 H.2 - cod. TAGMUA2 H.3 - cod. TAGMUA3 H.4 - cod. TAGMUA4</p>		<p>Angled Mua 17°</p> <p>H.2 - cod. TAGMUA172 H.3 - cod. TAGMUA173 H.4 - cod. TAGMUA174</p>		<p>Angled Mua 30°</p> <p>H.3 - cod. TAGMUA303 H.4 - cod. TAGMUA304 H.5 - cod. TAGMUA305</p>
	<p>Mua healing cap</p> <p>cod. TAGCG097</p>		<p>Mua transfer</p> <p>cod. TAGTRAMUA</p>		<p>Mua analog</p> <p>cod. TAGANAMUA</p>
	<p>Mua titanium temporary abutment</p> <p>cod. TAGMDMUA</p>		<p>Mua castable cylinder</p> <p>cod. TAGCCMUA</p>		<p>Abutment screw</p> <p>cod. TAGVM</p>
	<p>Mua screw M1.4</p> <p>cod. TAGVMMUA</p>		<p>Screwdriver</p> <p>L10 - cod. TAGDS L15 - cod. TAGDL</p>		<p>Ratchet screwdriver</p> <p>L12 - cod. TAGDCS L17 - cod. TAGDCL</p>
	<p>Machine screwdriver</p> <p>L20 - cod. TAGDMXS L26 - cod. TAGDMS L32 - cod. TAGDML</p>		<p>Ratchet implant driver</p> <p>L12 - cod. TAGIDCS L17 - cod. TAGIDCL</p>		<p>Machine implant driver</p> <p>Short - cod. TAGIDS Long - cod. TAGIDL</p>
	<p>Friction fit abutment extractor</p> <p>cod. TAGEXT</p>		<p>Mua machine screwdriver</p> <p>cod. TAGDMUA</p>		<p>Ratchet Mua screwdriver</p> <p>cod. TAGDCMUA</p>
	<p>Machine implant driver hand adapter</p> <p>cod. TAGPCM</p>		<p>Dynamometric torque ratchet adjustable 0 to 35 Ncm cod. TAGCRIDIN</p> <p>Fixed ratchet cod. TAGCRI</p>		<p>Cortical drill</p> <p>L26 - cod. TAGFL26 L32 - cod. TAGFL32</p>
	<p>DLC drill</p> <p>cod. TAGFF { 2.0 / 2.5 / 2.8 / 3.0 3.2 / 3.5 / 3.65 4.0 / 4.3 / 4.5 / 5.4</p>		<p>Drill extention</p> <p>cod. TAGAIP</p>		<p>Scan abutment</p> <p>cod. TAGSBI</p>

meté
Biomedical

Via Boccaccio, 8
21010 Arsago Seprio (VA)
Tel./Fax: +39 0331.796417
www.dentalmete.it